

# RX3 Vortex Mixer

## Vortex Mixer for Basic Requirements

- Fixed stirring speed, 3000 rpm
- Touch mode
- Enhanced stability on the bench
- New materials, higher performance



The RX3 is a basic vortex mixer, operating in **touch mode**.

**High-speed and excellent vortexing** is ensured by the fixed speed at 3000 rpm. The RX3 remains **extremely stable** thanks to special feet.

Get full speed by pressing down the rubber cup and mix your sample in an excellent way.

The new technopolymer/zinc alloy structure provides an **optimum chemical resistance**.

The RX3 is ideal for many applications, mixing samples in different containers.

VELP Scientifica offers many years of high performance in laboratories all over the world, with resistant and reliable solutions.

Live an amazing vortex experience!

## Features and Benefits

### Optimized Performance, Quiet Operation

Fixed speed (3000 rpm).

### Product Handling

Ergonomic and highly innovative design for excellent stability on the bench.

It is usable on many surfaces thanks to the 4 anti-sliding feet.

| Technical Data                  | Description                                              |
|---------------------------------|----------------------------------------------------------|
| Construction material:          | Zinc alloy and technopolymer                             |
| Support system:                 | 4 anti-sliding feet for strong fixing and high stability |
| Protection rating CEI EN 60529: | IP 42                                                    |
| Electronic speed regulation:    | constant, 3000 rpm                                       |
| Power Supply:                   | 100-240V/50-60Hz                                         |
| Power:                          | 15 W                                                     |
| Weight:                         | 2,7 Kg (5.9 lb)                                          |
| Dimensions (W x H x D) :        | 150x130x165 mm (5.9x5.1x6.5 in)                          |
| Ordering Information<br>Code No | Description                                              |
| F202A0171                       | RX3 Vortex Mixer 100-240V/50-60Hz                        |
| A00000016                       | Small rubber supporting plate Ø 50 mm                    |
| 10003083                        | US plug                                                  |
| 10003084                        | UK plug                                                  |
| 10003085                        | Australian plug                                          |



A00000016

Your authorized agent:



**tecno-lab** s.r.l.

Via L. Abbati, 22/A-B - 25131 Brescia - E-mail: [info@tecnolab.bs.it](mailto:info@tecnolab.bs.it)  
Tel. 0303582505 r.a. - Fax 0303582517 - [www.tecnolab.bs.it](http://www.tecnolab.bs.it)  
Apparecchiature scientifiche da laboratorio e assistenza tecnica

We reserve the right to make technical alternations  
We do not assume liability for errors in printing, typing or transmission