

VASCA COIBENTATA

Thermic insulation tank



Caratteristiche tecniche

- Vasche coibentate adatte a tutte le applicazioni di termostatazione, idonee per l'applicazione dei più comuni gruppi termostatici.
- Struttura esterna in acciaio trattata con vernice epossidica antiacido.
- Isolamento termico a doppia parete.
- Vasca interna in acciaio inox stampata con bordi arrotondati per agevolarne la pulizia.
- Rubinetto di scarico (a richiesta).
- Coperchio anti condensa in acciaio inox (a richiesta).
- Ponte per termostato in acciaio inox (a richiesta).
- Range di temperatura da - 20°C a +300°C (Serie VC).
- Range di temperatura da - 5°C a +120°C (Serie Basic).

Technical specifications

- Professional thermostatic tank ideal for the applications of most common thermostatic groups.
- Outer body in steel treated with anti-acid epoxy paint.
- Double wall heat insulation with natural mineral fibre.
- Internal chamber in seamless stainless steel with rounded corners for efficient circulation and cleaning.
- Draining tap (see accessories).
- Anti-condensation lid in stainless steel (see accessories).
- Bridge for thermostat stainless steel (see accessories).
- Temperature range from -20°C to +300°C (VC series).
- Temperature range from -5°C to +120°C (Basic series).

Modello / Model Codice / Code	M908 -VC 02.4001	M918 -VC 02.4002	M928 -VC 02.4003
Modello / Model Codice / Code	M908 -Basic 02.4101	M918 -Basic 02.4102	M928 -Basic 02.4103
DIMENSIONI INTERNE <i>INTERNAL DIMENSIONS</i>			
Volume nominale Lt. <i>Nominal volume Lt.</i>	8	18	28
Larghezza mm. <i>Width mm.</i>	300	330	503
Profondità mm. <i>Depth mm.</i>	150	298	298
Altezza mm. <i>Height mm.</i>	150	150	150
DIMENSIONI ESTERNE <i>OVERALL DIMENSIONS</i>			
Larghezza mm. <i>Width mm.</i>	325	360	535
Profondità mm. <i>Depth mm.</i>	175	325	325
Altezza mm. <i>Height mm.</i>	210	210	210
Porta provette inseribili max. (1) <i>Max insertable test tube rack</i>	0	1	2
Peso Kg. <i>Weight Kg.</i>	3	6	8



tecno-lab s.r.l.

Via L. Abbiati, 22/A-B - 25131 Brescia - E-mail: info@tecnolab.bs.it
Tel. 0303582505 r.a. - Fax 0303582517 - www.tecnolab.bs.it
Apparecchiature scientifiche da laboratorio e assistenza tecnica